**drive**

**DeVilbiss®**  
HEALTHCARE

**555 Series**  
**Serie de 555**  
**Série 555**

**en** DeVilbiss® 5-Liter Compact Oxygen Concentrator Instruction Guide

ASSEMBLED IN USA

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	<p><b>WARNING</b></p> <p>Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.</p>
--	---

**Physician Information:**

Physician Name: \_\_\_\_\_ Telephone: \_\_\_\_\_  
 Address: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_ Country: \_\_\_\_\_

**Prescription Information:**

Name: \_\_\_\_\_ Oxygen liters per minute \_\_\_\_\_ at rest: \_\_\_\_\_  
 during activity: \_\_\_\_\_ other: \_\_\_\_\_ Oxygen use per day Hours: \_\_\_\_\_ Minutes: \_\_\_\_\_  
 Comments: \_\_\_\_\_

**DeVilbiss 5-Liter Oxygen Concentrator Serial Number:** \_\_\_\_\_

**DeVilbiss Equipment Provider Information:**

Set-Up Person: \_\_\_\_\_

This instruction guide was reviewed with me, and I have been instructed on the safe use and care of the DeVilbiss Oxygen Concentrator

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

# WARNINGS AND PRECAUTIONS

## IMPORTANT SAFEGUARDS

Read this entire guide before using your DeVilbiss concentrator. Important safeguards are indicated throughout this guide. Pay special attention to all safety information. Imminently and potentially hazardous information is highlighted by these terms:

### DANGER

Indicates an imminently hazardous situation which could result in death or serious injury to the user or operator if not avoided.

### WARNING

Indicates a potentially hazardous situation which could result in death or serious injury to the user or operator if not avoided.

### CAUTION

Indicates a potentially hazardous situation which could result in property damage, injury, or device damage if not avoided.

### IMPORTANT

Indicates important information you should know.

### NOTE

Indicates notes, useful tips, recommendations, and information.

## READ ALL INSTRUCTIONS BEFORE USE.

### DANGER

- NO SMOKING signs should be prominently displayed when the concentrator is in use.
- Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy.
- Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located.
- If you intend to smoke, you must always turn the oxygen concentrator OFF, remove the cannula and leave the room where either the cannula or mask or the oxygen concentrator is located. If unable to leave the room, you must wait 10 minutes after you have turned OFF the oxygen concentrator before smoking.
- Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned ON but not in use. The oxygen will make the materials flammable. Turn the oxygen concentrator OFF when not in use to prevent oxygen enrichment.
- Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.
- Open flames during oxygen therapy are dangerous and are likely to result in fire or death. Do not allow open flames within 2 m (6.5 feet) of the oxygen concentrator or any oxygen carrying accessories.
- DeVilbiss oxygen concentrators are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit.

### WARNING

- To avoid electric shock do not plug the concentrator into an AC outlet if the concentrator exterior cabinet is broken. Do not remove the concentrator exterior cabinet. The exterior cabinet should only be removed by a qualified service technician. Do not apply liquid directly to the exterior cabinet or utilize any petroleum-based solvents or cleaning agents.
- Improper use of the power cord and plugs can cause a burn fire or other electric shock hazards. Do not use the unit if the power cord is damaged.
- Ensure the mains power cord is fully inserted into the concentrator connector and the power cord plug is completely inserted into a fully functioning AC wall outlet. Failure to do so, may cause an electrical safety hazard.
- The accessories (nasal cannula, masks, oxygen tubing, humidifiers) that supply oxygen to the patient must be equipped with a means that in case of fire, stops the propagation of fire through the accessory for the safety of the patient and others. A fire activated flow-stop or thermal fuse device if available should be used with the oxygen supply accessories. These types of flow-stop devices stop the flow of oxygen to the patient in the event of fire. This means of fire protection should be located as close to the patient as practicable.
- Locate oxygen tubing and power supply cords to prevent tripping hazards and reduce the possibility of entanglement or strangulation.
- Do not lubricate fittings, connections, tubing or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
- Do NOT use lubricants, oils or grease.
- Before attempting any cleaning procedures, turn the unit "OFF".
- Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.
- Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
- Keep the oxygen concentrator out of the reach of children and pets to avoid a change or loss of therapy.
- This device may be affected by children unintentionally changing the device setting or inserting foreign objects into device openings. Supervise children when the device is in operation.
- The device can be affected by pet hair or insects getting into the device or obstructing cooling vents in the device.

## WARNING

- If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy seek medical assistance immediately to avoid harm.
- Geriatric, pediatric or any other patient unable to communicate discomfort can require additional monitoring and/or a distributed alarm system to convey the information about the discomfort and/or the medical urgency to the responsible caregiver to avoid harm.
- Use of this device at an altitude above 9840 feet (3000 meters) or above a temperature of 95°F (35°C) or greater than 90% relative humidity may affect the flow rate and the percentage of oxygen and consequently the quality of the therapy. Refer to specifications for details regarding parameters tested.
- To ensure you receive the therapeutic amount of oxygen delivery according to your medical condition the Oxygen Concentrator must:
  - be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
  - be used with the specific combination of parts and accessories that are in line with the specification of the concentrator manufacturer and that were used while your settings were determined.
- Your delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.
- For your safety the oxygen concentrator must be used according to the prescription determined by your physician.
- Under certain circumstances oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.

## WARNING

### MR UNSAFE

- Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the oxygen concentrator or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.
- Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the oxygen concentrator. Some electromagnetic sources may not be apparent. If you notice any unexplained changes in the performance of this device or if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your home care provider.
- This device is suitable for use in home and healthcare environments except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging where the intensity of Electromagnetic DISTURBANCES is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals, such as antenna cables and external antennas) should be used no closer than 30.5 cm (12 inches) to any part of the oxygen concentrator including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

## WARNING

### Risk of injury or damage

- When the device is operated at the extremes of the environmental operating specifications (i.e. maximum temperatures and humidity) and in a single fault condition which is a single component or performance malfunction such as a blocked exhaust vent blocked air intake or an internal cooling fan failure the temperature of the air coming out of the exhaust vents located on the right side of the unit can reach temperatures capable of causing a burn injury (see stated temperature and contact time values in the model specific tables below).
- Keep exposed body parts such as hands and feet a minimum of 12 inches (30.5 cm) away from the exhaust vents to avoid the risk of burns. Single fault conditions may result in visual and audible alerts and alarms.
- **NOTE** – Under normal and single fault conditions the concentrator releases warm air out the bottom of the unit (exhaust vents) which may discolor temperature sensitive flooring surfaces. The concentrator should not be used over flooring that is sensitive to heat staining. The Manufacturer is not responsible for flooring that becomes discolored.
- Operate the unit in a cool dry area with good ventilation located on a hard surface, avoid thick rugs or carpeting. NEVER block the air intake or exhaust vents. Keep the unit a minimum of 12 inches (30.5 cm) away from any wall draperies or any other objects that might prevent the proper flow of air in and out of your oxygen concentrator. Proper air flow is needed to prevent overheating of the oxygen concentrator. DO NOT place the concentrator near any heat source, such as hot air registers or heaters. Overheating of the oxygen concentrator may lead to low oxygen output and a risk of burns.
- The air intake and the exhaust of the oxygen concentrator should be located in a well-ventilated area. DO NOT operate the unit in a closed or confined space such as a closet, bathroom, etc. Avoid operating the device near smoke pollutants and fumes.
- If under extreme environmental conditions and a single fault condition occurs, the following device surface temperatures may exceed 106 °F (41 °C). See Table 1 below for the model specific maximum temperature and safe contact guidance:

Table 1 - 555DS Series

Description	Maximum Temperature		Max safe contact time
	555DS Series		
	°F	°C	
Oxygen Outlet Fitting	140.0	60.0	Less than 1 minute
Cannula at the Oxygen Outlet Fitting	112.8	44.9	Less than 10 minutes
Power Switch	121.8	49.9	Less than 10 minutes
Exhaust Vent AIR	183.6	84.2	Less than 1 second
Flow Control Knob	121.5	49.7	Less than 10 minutes
Exhaust Vent Plastic Fin	153.0	67.2	Less than 10 seconds
Exterior Cabinet (rear cover)	128.0	53.3	Less than 10 minutes
Circuit Breaker External Surface	133.3	56.3	Less than 1 minute
IEC Power Connector	127.0	52.8	Less than 10 minutes

See Important Parts of Your Concentrator in this instruction guide for location of accessible parts listed in the table.

## CAUTION

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.
- Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage.

## CAUTION

- To avoid possible injury or loss of therapy, turn the device off and disconnect from patient before moving the device.



## IMPORTANT

- Do not service or clean this device while in use with a Patient.
- The Device is classified as IP21 which means it is protected against finger access to hazardous parts and protected against vertically falling water drops.
- This device contains electrical and/or electronic equipment. Follow local governing ordinances and recycling plans regarding disposal of device components.
- Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

# SAVE THESE INSTRUCTIONS

## SYMBOLS GLOSSARY

 <p>It is mandatory to read and understand the operating instructions prior to use.</p> <p><b>i</b> This symbol has a blue background on the product label.</p>	 <p>OFF ON</p>	 <p>Operating Temperature Range +5 to +35°C (+41 to +95°F)</p>	 <p>Atmospheric Pressure Range 715 to 1013 hPa (Approximate sea level to 9840 ft)</p>
 <p>Electric Shock Hazard. Exterior Cabinet to be removed by authorized personnel only.</p> <p><b>i</b> This symbol has a yellow background on the product label.</p>	 <p>Manufacturer</p>	 <p>Transport and Storage Relative Humidity Range 15% to 90%, non-condensing</p>	 <p>Transport and Storage Temperature Range -13°F (-25°C) to 158°F (70°C)</p>
 <p>Danger - No smoking near patient or device.</p> <p><b>i</b> This symbol has a red circle and diagonal bar on the product label.</p>	 <p>Alternating Current</p>	 <p>Serial Number</p>	 <p>European CE mark</p>
 <p>Use no Oil, Grease or Lubricants</p> <p><b>i</b> This symbol has a red circle and diagonal bar on the product label.</p>	 <p>Type BF applied part</p>	 <p>Medical Device</p>	 <p>Keep unit dry.</p>
 <p>Do not use near heat or open flames</p> <p><b>i</b> This symbol has a red circle and diagonal bar on the product label.</p>	 <p>Class II Equipment</p>	 <p>Normal Oxygen</p> <p><b>i</b> This symbol is green on a black background.</p>	 <p>Ingress Protection - Protected against finger access to hazardous parts; protected against vertically falling water drops.</p>
 <p>Do Not Disassemble</p> <p><b>i</b> This symbol has a red circle and diagonal bar on the product label.</p>	 <p>Note and Information Symbol</p> <p><b>i</b> This symbol is used throughout this manual to indicate notes, useful tips, recommendations and information.</p>	 <p>Oxygen Outlet marked with rated range of gas pressure and rated range of gas flow rate.</p>	 <p>This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU- Waste Electrical and Electronic Equipment. (WEEE)</p>
 <p>MRI Unsafe Resonance Environment This symbol has a red circle and diagonal bar on the product label.</p>	 <p>Catalog Number</p>	 <p>Low Oxygen -</p> <p><b>i</b> This symbol is yellow on a black background.</p>	 <p>CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.</p>
 <p>General Warning</p> <p><b>i</b> This symbol is used throughout this manual to indicate hazardous situations to avoid.</p>	 <p>European Authorized Representative</p>	 <p>Consult Instructions for use</p>	 <p>Unique Device Identifier</p>
 <p>Important Information</p> <p><b>i</b> This symbol is used throughout this manual to indicate important information you should know.</p>	 <p>LOT Number</p>	 <p>Service Required</p> <p><b>i</b> This symbol is red with a black background.</p>	 <p>Flow Rate Display Maximum recommended flow rate: 5LPM</p>
 <p>Startup Period Low O2 (&lt; 85.0 %)</p> <p><b>i</b> This symbol is yellow on a black background.</p>	 <p>Hour meter</p>	 <p>TUV Rheinland C-US approval mark.</p>	 <p>Warning – O2S gas over temperature.</p> <p><b>i</b> This symbol is red on a black background.</p>
 <p>Swiss Authorized Representative</p>	 <p>Caution is necessary when operating the device.</p>	 <p>Increase flow clockwise Decrease flow counter-clockwise.</p>	

## **INTRODUCTION**

This instruction guide will provide you with information on the use of your DeVilbiss Oxygen Concentrator. Make sure that you read and understand this guide before operating your unit. Important safeguards are indicated throughout this guide. Pay special attention to all safety information. Contact your DeVilbiss equipment provider should you have any questions.

## **INTENDED USE**

The DeVilbiss 555 5L Compact Oxygen Concentrator is intended to provide supplemental oxygen to adult patients requiring oxygen therapy.

## **CONTRAINDICATIONS**

None known.

## **INDICATIONS FOR USE**

The DeVilbiss 555 5L Compact Oxygen Concentrator is intended to provide supplemental oxygen to adult patients requiring oxygen therapy. The device may be used in the home or an institutional setting. The device is not intended to be life-sustaining or life-supporting.

## **ESSENTIAL PERFORMANCE**

Essential Performance of the Oxygen Concentrator is to deliver a continuous flow of oxygen enriched gas. Visual and audible alarms indicate if the device is not meeting specification or a failure has been detected. If Essential Performance of the oxygen concentrator is lost or degraded, the device is producing less than 85% oxygen, resulting in a low O<sub>2</sub> alarm.

## **SERVICE LIFE**

The expected service life of the 555 series 5L Compact Oxygen Concentrator, which includes the performance of any required service or maintenance, is 5 years. The expected service life is based on the operation of the device in accordance with all manufacturer guidance for safe use, maintenance, servicing, storage, shipping, handling, and general operation.

The actual service life of the unit, and in particular the service life of certain subcomponents, including the Filters, Sieve Beds and Compressor Cup Seals, will vary based on a number of variables, including the operating environment, storage environment, shipping, handling, performance of preventive maintenance, and both the frequency and intensity of use.

The 555 series 5L Compact Oxygen Concentrator have internal sensors and diagnostic systems designed to monitor the system performance, including the oxygen concentration (purity), flow and temperature. The 555 series 5L compact oxygen concentrator will alert the user when the device requires maintenance or service. Please see the Troubleshooting and Maintenance Sections for more detailed information.

### **NOTE**

Drive DeVilbiss recommends for optimal service life that the 555 Series 5L Compact Oxygen Concentrator be operated for at least 30 minutes after it is powered on. Shorter periods of operation, operating in extreme temperature/humidity conditions or in the presence of contaminants, and/or handling and storage conditions outside those specified, may affect the long term reliable operation of the product.

## **WHY YOUR PHYSICIAN PRESCRIBED SUPPLEMENTAL OXYGEN**

Today, many people suffer from heart, lung and other respiratory diseases. Many of these people can benefit from supplemental oxygen therapy. Your body requires a steady supply of oxygen to function properly. Your physician prescribed supplemental oxygen for you, because you are not getting enough oxygen from room air alone. Supplemental oxygen will increase the amount of oxygen that your body receives.

Supplemental oxygen is not addictive. Your physician prescribed a specific oxygen flow to improve symptoms such as headaches, drowsiness, confusion, fatigue or increased irritability. If these symptoms persist after you begin your supplemental oxygen program, consult your physician.

The oxygen delivery setting has to be determined for each patient individually with the configuration of the equipment to be used, including accessories.

The proper placement and positioning of the prongs of the nasal cannula in the nose is critical to the amount of oxygen delivered to the respiratory system of the patient.

Your Delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.

## **HOW YOUR DeVilbiss OXYGEN CONCENTRATOR WORKS**

Oxygen concentrators are the most reliable, efficient and convenient source of supplemental oxygen available today. The oxygen concentrator is electrically operated. The unit separates oxygen from room air which allows high-purity supplemental oxygen to be delivered to you through the oxygen outlet. Although the concentrator filters the oxygen in a room, it will not affect the normal amount of oxygen in your room. If the oxygen purity falls below the acceptable level, the display will show a Yellow Low Oxygen symbol, and an intermittent audible signal will sound.

If the oxygen purity continues to fall into the critical low level, the Red Wrench (Service Required) symbol will show on the display and an intermittent audible signal will sound. Refer to the Troubleshooting section in this guide. Do not attempt any other maintenance. Contact your service technician immediately if the issue is not resolved.

### **IMPORTANT**

- Do not service or clean this device while in use.
- The 555 Series 5L Compact Oxygen Concentrator, its parts and accessories are specified for use at flow rates between 0.5 LPM and 5 LPM.

## **ACCESSORIES**

Bubble Humidifier .....Salter Labs 7600 or equivalent

High Flow Nasal Cannula.....Salter Labs 1600HF or equivalent

Recommended Filter

Compressor intake filter.....Drive DeVilbiss 555D-615

Firesafe Cannula Valve.....RES010

Refer to the accessory manufacturer's instructions for use. A mask or any nasal cannula can be used with continuous flow delivery and may be sized according to your prescription as recommended by the clinical staff who should also give you advice on the proper usage, maintenance, and cleaning.

### **WARNING**

The accessories (nasal cannula, masks, oxygen tubing, humidifiers, etc.) that supply oxygen to the patient must be equipped with a means that, in case of fire, stops the propagation of fire through the accessory for the safety of the patient and others. A fire activated flow-stop or thermal fuse device, if available, should be used with the oxygen supply accessories. These types of flow-stop devices stop the flow of oxygen to the patient in the event of fire. This means of fire protection should be located as close to the patient as practicable.

- ▶ Do not use accessories, detachable parts and materials not described in the instructions for use.
- ▶ Do not modify the oxygen concentrator or accessories.
- ▶ It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.
- ▶ There are many types of humidifiers, oxygen tubing and cannulas/masks that can be used with this device. Certain humidifiers and accessories may impair the device's performance.
- ▶ Do not connect this oxygen concentrator with other equipment not described in the instructions for use.

### **CAUTION**

- Incompatible parts or accessories can result in degraded performance.
- Humidifier accessories may create a hazard if they are overfilled or overflow onto the oxygen concentrator. Excessive liquid poured over the oxygen concentrator may cause the oxygen concentrator to malfunction.


### **NOTE**

- The device should only be used with bubble humidifiers that are designed for use with flows up to 5 liters per minute and 8.5 psi pressure.
- A maximum of 50 feet (15.2 meters) of crush-proof oxygen tubing, plus 7 feet (2.1 meters) of cannula, plus a bubble humidifier is allowed between the concentrator and the patient.
- The clinical staff or your provider should verify the compatibility of the oxygen concentrator and all the parts used to connect to the patient before use.
- The clinical staff or your provider may have preset the oxygen flow so that it cannot be adjusted.
- Your physician has prescribed either a nasal cannula or face mask. In most cases, they are already attached to the oxygen tubing. If not, follow the manufacturer's instructions for attachment.

## IMPORTANT PARTS OF YOUR CONCENTRATOR

Please familiarize yourself with your DeVilbiss Oxygen Concentrator before use.

### FRONT/SIDE VIEW (FIGURE A)

1. Oxygen outlet fitting
2. Control Panel
3. Handle
4. Exhaust Vents
5. Casters

### BACK VIEW (FIGURE B)

6. Humidifier Bottle Strap
7. Humidifier Bottle Holder
8. Line cord strap
9. Intake Filter Door
10. Intake Filter (inside of Intake Filter Door)
11. IEC power connector
12. Circuit Breaker - resets the unit after electrical overload shutdown

### CONTROL PANEL DETAIL (FIGURE C)

13. Digital Display
14. Flow adjustment knob
15. Power Switch

| = ON

O = OFF

### WARNING

When the device is used under extreme operating or single fault conditions the exhaust air temperature near the exhaust vents on the bottom of the unit may exceed 106°F (41°C). Keep exposed body parts such as hands and feet a minimum of 12 inches (30.5 cm) away from the exhaust vents to avoid the risk of burns.

### Location of the Intake Filter



### To Attach/Detach V-Lock Power Cord

1. To attach the power cord, insert the power cord into the IEC connector on the concentrator. The end of the power cord clicks when it locks in place.

2. To detach the power cord, a) press the yellow tab on the power cord connector. b) Pull the power cord from the IEC connector on the concentrator.

**NOTE** - The Power cord locks in place. Do not forcefully try to remove the power cord; the locking tab must first be pressed in to remove the power cord.


## **SETTING UP YOUR OXYGEN CONCENTRATOR**

1. With the power switch of the device in the OFF position (Figure A), insert the power cord into the IEC Connector located on the back side of the device (Figure B). Insert the other end of the power cord into the wall outlet as shown below (Figure C).

**A**

**B**


**C**


2. Position your unit near an electrical outlet in the room where you spend most of your time.
3. Position your unit on a flat surface at least 12 inches (30.5 cm) from walls, draperies or any other objects that might prevent the proper flow of air in and out of your oxygen concentrator. The oxygen concentrator should be located in a well-ventilated area to avoid pollutants or fumes.
4. Before operating your unit, always check to be sure the filter door and exhaust vents (located on the sides of your unit) are clean. Proper cleaning is discussed in the Caring For Your DeVilbiss Oxygen Concentrator section of these instructions for use.
5. Attach the oxygen tubing to the oxygen outlet connector. (Figure D)

**D**

## **OXYGEN TUBING CONNECTION WITH HUMIDIFICATION**

1. Follow steps 1 through 4 above to connect the device to AC power and position it for use.
2. If your physician has prescribed an oxygen humidifier as part of your therapy, follow these steps (If using a prefill, go to step b.):
  - a. Fill the humidifier bottle with distilled water.
  - b. Insert the humidifier into the humidifier holder and secure the humidifier in place using the humidifier strap. (Figure E)
  - c. Thread the wing nut located on the top of the humidifier bottle to the humidifier tubing. Make sure it is securely tightened. (Figure F)
  - d. Attach the other end of the humidifier tubing to the oxygen outlet fitting on the 555 Series 5L Compact Oxygen Concentrator. (Figure G)
  - e. Attach the patient oxygen tubing directly to the humidifier bottle outlet fitting. (Figure H)

	<b>WARNING</b> Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged.
--	--

	<b>WARNING</b> Ensure the mains power cord is fully inserted into the concentrator connector and the power cord plug is completely inserted into a fully functioning AC wall outlet. Failure to do so may cause an electrical safety hazard.
--	--

- i** **NOTE** – Do not connect to an electrical outlet controlled by a wall switch.
- i** **NOTE** – To move the unit, firmly grasp the handle located on the top of the unit, rolling and/or lifting the unit over pathway obstacles.
- i** **NOTE** – Your physician has prescribed either a nasal cannula or face mask. In most cases, they are already attached to the oxygen tubing. If not, follow the manufacturer's instructions for attachment.
- i** **NOTE** – Your healthcare provider should verify the compatibility of the oxygen concentrator and all the parts used to connect to the patient before use.
- i** **NOTE** – (only 120 volt units) The plug on the DeVilbiss oxygen concentrator has one blade wider than the other. To reduce the risk of electric shock, this plug is intended to fit in a wall outlet only one way. Do not attempt to defeat this safety feature.
- i** **NOTE** – To check your oxygen concentrator and accessories for proper operation; 1. Check the output flow by placing the end of the nasal cannula under the surface of a half-full cup of water and look for bubbles. 2. Check the system for leaks by bending the nasal prongs over and squeeze tight to stop the flow of oxygen. Look at the LCD display to see that the flow rate on the display drops to zero. If the displayed flow rate does not drop to zero, check all connections for possible leaks. Parts to check for leaks are: tubing connections, humidifier bottle and other accessories like firebreaks. Repeat these steps until the displayed flow rate drops to zero. 3. Contact your provider or service supplier immediately if you encounter any problems.

## **OPERATING YOUR DeVilbiss OXYGEN CONCENTRATOR**

1. Press the power switch to the “ON” position. When the unit is turned “ON,” the LCD display shows a Drive DeVilbiss logo and the hour meter. An audible signal will briefly alarm confirming that the LCD display and audible signal are functioning properly. The unit operates in “start up” mode with the Startup Period Low O<sub>2</sub> symbol shown on the LCD display until a normal oxygen level is achieved. When a normal oxygen level is achieved, the Startup Period Low O<sub>2</sub> symbol turns off and the Normal Oxygen symbol shows on the LCD display. The “start up” period may take up to 15 minutes.
2. Check the flow setting shown on the LCD display to make sure it matches your prescribed oxygen flow rate.

3. Example: LCD display showing 5.0 Liters per minute flow setting

4. Your DeVilbiss concentrator is now ready for use, properly position the cannula with the nasal prongs facing upward, insert the prongs into nose. Wrap the cannula tubing over the ears and position in front of body (Figure I).

	<p><b>DANGER</b></p> <ul style="list-style-type: none"><li>• Oxygen causes rapid burning.</li><li>• Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy.</li><li>• Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.</li><li>• DeVilbiss oxygen concentrators are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit.</li></ul>	<p><b>NO SMOKING</b></p>
--	---	------------------------------

	<p><b>DANGER</b></p> <ul style="list-style-type: none"><li>• Oxygen makes it easier for a fire to start and spread.</li><li>• Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned ON but not in use, the oxygen will make the materials flammable.</li><li>• Turn the oxygen concentrator OFF when not in use to prevent oxygen enrichment.</li></ul>
--	--

	<p><b>WARNING</b></p> <p>In order to prevent a fire propagating from the patient through the cannula towards the unit, a means of protection should be located as close to the patient as practical. Please contact your provider for this means of protection.</p>
--	---

	<p><b>CAUTION</b></p> <p>When the unit is turned “ON”, as part of the normal start-up process, the LCD Display should appear on the digital display and the audible alarm should sound. If the digital display DOES NOT turn on or the audible alarm DOES NOT sound, this indicates the alert system is not functioning properly. Refer to the Troubleshooting chart and contact your DeVilbiss oxygen provider if necessary.</p>
--	---

	<p><b>CAUTION</b></p> <p>It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.</p>
--	--

**i NOTE**– DeVilbiss recommends for optimal service life that the DeVilbiss Oxygen Concentrator be operated for at least 30 minutes after it is powered ON. Shorter periods of operation, operating in extreme temperature/humidity conditions or in the presence of contaminants, and/or handling and storage conditions outside those specified, may affect the long term reliable operation of the product.

**i NOTE** – If the audible signal alarms but the unit is not operating, there is no power to the unit. Refer to the Troubleshooting chart and contact your DeVilbiss provider if necessary.

- i** **NOTE** – If an audible low-frequency vibration sound is detected, the unit is not operating properly. Refer to the Troubleshooting chart and contact your DeVilbiss provider if necessary.
- i** **NOTE** – If the flow adjustment knob is turned counterclockwise, the flow decreases (and eventually will shut off the oxygen flow). If the knob is turned clockwise, the flow increases.
- i** **NOTE** – For prescriptions of 5 LPM, be sure the display shows 5.0. Setting the flow higher than 5.0 may cause the oxygen purity level to drop, and Yellow O2 with Down Arrow symbol will show on the LCD display, accompanied by an audible alarm. Adjust the oxygen flow to your prescribed flow rate.
- i** **NOTE** – The low-flow alarm may activate if the oxygen flow is set at or below 0.2 LPM. The unit will continue to run; however, the Red Wrench Service Required symbol will show on the LCD display, accompanied by an audible alarm. Adjust the oxygen flow to your prescribed flow rate.

## **DeVilbiss OSD® OPERATION**

The OSD (Oxygen Sensing Device) is a device within your concentrator that monitors the oxygen produced by your unit. The oxygen status symbols on the LCD display are defined as follows:

Green O2 with check mark symbol on the LCD display indicates output is an acceptable oxygen level.

Yellow O2 with Down Arrow symbol on the LCD display indicates output is below an acceptable oxygen level.

If the oxygen purity falls below the acceptable level: The Green O2 with check mark symbol will shut off, the Yellow O2 with Down Arrow symbol will illuminate, and an intermittent audible signal will sound.

Refer to the Troubleshooting section in this guide and switch to your reserve oxygen system. Do not attempt any other maintenance. Contact your DeVilbiss provider immediately.

## **RESERVE OXYGEN SYSTEM**

As a precaution, your DeVilbiss provider may supply you with a reserve oxygen system. If your oxygen concentrator loses electrical power or fails to operate correctly, the oxygen concentrator's alarm system will sound to signal you to switch to your reserve oxygen system (if provided) and contact your DeVilbiss provider. Please contact your oxygen provider if you have questions regarding a reserve oxygen system.

### **DANGER**

Before attempting any cleaning procedures, turn the unit OFF and disconnect from mains power. Risk of electrical shock.

## **CARING FOR YOUR DeVilbiss OXYGEN CONCENTRATOR**

DeVilbiss recommends using only original DeVilbiss parts and filters to guarantee reliable operation of the product.

### **CANNULA/MASK, TUBING AND HUMIDIFIER BOTTLE**

Thoroughly clean or replace the cannula/mask, tubing and humidifier bottle according to the manufacturer's instructions.

### **FILTER DOOR WITH VENTS**

Inspect the filter door vents periodically, and wipe with a dry cloth to remove dust.

### **EXHAUST VENTS**

Inspect the exhaust vents periodically, and wipe with a dry cloth to remove dust.

### **EXTERIOR CABINET**

Thoroughly clean the concentrator exterior cabinet weekly or if the exterior cabinet is dusty or dirty by using a damp cloth and wiping dry; the vents can also be wiped with a damp cloth.

### **CLEANING INSPECTION**

After cleaning, visually inspect the parts. If they are not visually clean, repeat the previous cleaning steps or dispose of so that a visibly soiled device is not reused. Allow cleaned parts to dry completely before using again.

## CLEANING

	<b>RECOMMENDED CLEANING INTERVAL</b>	<b>NUMBER OF CLEANING CYCLES *</b>	<b>COMPATIBLE CLEANING METHOD</b>
Exterior Cabinet	7 days	260	Clean tap water, use only a damp cloth
Filter Door and Exhaust Vents	7 days	260	Wipe with dry cloth, or a cloth dampened with clean tap water to remove dust
Oxygen Outlet Connector	7 days	104	Mild dish soap (2 tbsp) and clean warm tap water (2 cups)

\* number of cleaning cycles determined by recommended cleaning interval and expected service life

	<p><b>WARNING</b></p> <ul style="list-style-type: none"> <li>Do NOT use lubricants, oils or grease.</li> <li>Before attempting any cleaning procedures, turn the unit “OFF.”</li> </ul>
--	---

	<p><b>CAUTION</b></p> <p>Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage.</p>
--	---

## TROUBLESHOOTING

The following troubleshooting chart will help you analyze and correct minor oxygen concentrator malfunctions. If the suggested procedures do not help, switch to your reserve oxygen system and call your DeVilbiss home care provider. Do not attempt any other maintenance.

	<p><b>WARNING</b></p> <ul style="list-style-type: none"> <li>To avoid electric shock, do not plug the concentrator into an AC outlet if the concentrator exterior cabinet is broken. Do not remove the concentrator exterior cabinet.</li> <li>The exterior cabinet should only be removed by a qualified DeVilbiss technician.</li> </ul>
---	--

## **TROUBLESHOOTING CHART**

<b>SYMPTOM</b>	<b>POSSIBLE CAUSE</b>	<b>REMEDY</b>	<b>ACTION IF REMEDY FAILS TO RESOLVE ISSUE</b>
<p><b>A. Unit does not operate.</b></p> <p>Display is off when the power switch is “ON.”</p> <p>Audible alert is pulsing.</p>	1. Power cord not properly inserted into wall outlet.	1. Check power cord connection at the wall outlet, also check the power cord connection on the back of the unit.	<p><b>If NONE of these remedies fix the issue, contact your Home Oxygen Provider, the device may require servicing.</b></p>
	2. No power at wall outlet.	2. Check your home circuit breaker and reset if necessary, Use a different wall outlet if the situation occurs again.	
	3. Oxygen concentrator circuit breaker activated.	3. Press the concentrator circuit breaker reset button located above the power cord connection. Use a different wall outlet if the situation occurs again.	
<p><b>B. Unit operates.</b></p> <p>Red Service Required symbol is illuminated.</p> <p>Audible alert may be sounding.</p>	1. Filter door vents are blocked.	1. Check filter door vents and ensure that the openings are not blocked.	<p><b>If NONE of these remedies fix the issue, contact your Home Oxygen Provider, the device may require servicing.</b></p>
	2. Exhaust is blocked.	2. Check the exhaust area and make sure there is nothing restricting the unit exhaust.	
	3. Blocked or defective cannula, face mask, or oxygen tubing.	3. Detach cannula or face mask. If proper flow is restored, clean or replace, if necessary. Disconnect the oxygen tubing at the oxygen outlet. If proper flow is restored, check oxygen tubing for obstructions or kinks. Replace if necessary.	
	4. Blocked or defective humidifier bottle.	4. Detach the humidifier from the oxygen outlet. If proper flow is obtained, clean or replace humidifier.	
	5. Flow meter set too low.	5. Set flow meter to prescribed flow rate.	

SYMPTOM	POSSIBLE CAUSE	REMEDY	ACTON IF REMEDY FAILS TO RESOLVE ISSUE
C. Unit operates. Audible low-frequency vibration sound is detected.	1. Electronic Assembly Malfunction.	1. Turn your unit "OFF " Switch to your reserve oxygen system and contact your DeVilbiss provider immediately.	If this remedy DOES NOT fix the issue, contact your Home Oxygen Provider, the device may require servicing.
D. The yellow Low Oxygen symbol is on.	1. Unit is in "start up" mode.	1. Allow unit up to 15 minutes to complete start up period.	If this remedy DOES NOT fix the issue, contact your Home Oxygen Provider, the device may require servicing.
E. The yellow Low Oxygen symbol is on, and the intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number.	If NONE of these remedies fix the issue, contact your Home Oxygen Provider, the device may require servicing.
	2. Filter door vents are blocked.	2. Check filter door vents and ensure that the openings are not blocked.	
	3. Exhaust is blocked.	3. Check the exhaust area and make sure there is nothing restricting the unit exhaust.	
F. Red Service Required symbol is on and an intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number.	If NONE of these remedies fix the issue, contact your Home Oxygen Provider, the device may require servicing.
	2. Filter door vents are blocked.	2. Check filter door vents and ensure that the openings are not blocked.	
	3. Exhaust is blocked.	3. Check the area near the exhaust vents and make sure there is nothing restricting the unit exhaust.	
	4. Electronic Assembly Malfunction.	4. Turn your unit "OFF " Switch to your reserve oxygen system and contact your DeVilbiss provider immediately.	
G. If any other problems occur with your oxygen concentrator.		1. Turn your unit "OFF." Switch to your reserve oxygen system and contact your oxygen provider immediately.	Contact your Home Oxygen Supplier, the unit requires servicing.
H. Unit operates. Any of the visual and audible alerts do not function when the power switch is turned "ON."	1. Electronic assembly malfunction.	1. Turn your unit "OFF." Switch to your reserve oxygen system and contact your oxygen provider immediately.	Contact your Home Oxygen Supplier, the unit requires servicing.
I. Red Service Required Symbol is on and there is an intermittent audible signal sounding. Unit may shut down.	1. Unit may be operating above normal temperature range.	1. Turn your unit "OFF." Switch to your reserve oxygen system and contact your oxygen provider immediately.	If NONE of these remedies fix the issue, contact your Home Oxygen Supplier, the unit may require servicing.
		2. Ensure the unit is not operating too close to a heat source, such as a fireplace, furnace duct/outlet, portable heater.	
		3. Ensure the unit is away from walls, curtains, etc., and is not in an enclosed space, such as a very small room or closet.	

## OVERVIEW OF ALARMS AND SERVICE INDICATORS

This device contains an alarm system which monitors the state of the device and alerts of abnormal operation, loss of essential performance or failures. Alarm conditions are shown on the LCD display. The alarm system functions are tested at power up by lighting all visual alarm indicators and sounding the audible alarm (beep).

All alarms are Low Priority Technical Alarms.

ALERT OR ALARM CONDITION	LCD ICON	DETAILS OF ALERT OR ALARM CONDITION	VISUAL ALERT OR ALARM	AUDIBLE ALARM	ACTION
Start-up Period  This symbol is yellow with a black background		The unit has recently been started and is in start-up period, the output flow of the oxygen is temporarily < 85%.	The Startup Period Low O2 symbol on the LCD is illuminated indicating low O2 condition.	No audible alarm during start-up period.	Wait for unit to finish start-up period, up to 15 minutes.

ALERT OR ALARM CONDITION	LCD ICON	DETAILS OF ALERT OR ALARM CONDITION	VISUAL ALERT OR ALARM	AUDIBLE ALARM	ACTION
Low Oxygen Output Concentration  This symbol is yellow with a black background		The output flow of oxygen is < 85%, which indicates the unit may need routine servicing.	The Low O2 symbol on the LCD is illuminated, indicating a Low O2 condition.	The audible alarm is beeping intermittently.	Contact your Oxygen Equipment Provider for assistance and to arrange for servicing of the unit.
High Temperature  This symbol is red with a black background		Internal device temperature is over the safe limit.	The RED Wrench Service Required with Thermometer symbol on the LCD is illuminated.	The audible alarm is beeping intermittently.	Turn the concentrator OFF, allow the device to cool. If the problem persists, contact your Oxygen Equipment Provider for assistance and to arrange for servicing of the unit.
Device Malfunction  This symbol is red with a black background		The device is experiencing a malfunction that requires servicing to correct. These malfunctions may include O2 < 60%, high flow, obstructed flow, fan not running, nonrecoverable valve error, corrupted settings, O2S sensor error.	The RED Wrench Service Required symbol on the LCD is illuminated.	The audible alarm is beeping intermittently.	Contact your Oxygen Equipment Provider for assistance and to arrange for servicing of the unit.

## **SPECIFICATIONS**

CATALOG NUMBER	555DS
Delivery Rate	0.5 to 5 LPM
Maximum Recommended Flow (@ nominal outlet pressures of zero & 7 kPa)	5 LPM
Outlet Pressure	1 - 5 LPM: 8.5 ± 0.5 PSIG (58.6 ± 3.5 kPa) 0.5 - < 1 LPM: 8.5 ± 0.75 PSIG (58.6 ± 5.2 kPa)
Electrical Rating	120 V~, 60 Hz, 2.6 A Max
Operating Voltage Range	102-132 V~, 60 Hz
Oxygen Percentage	0.5-5 LPM = 87%-96% O2
Operating Temperature Range	41°F (5°C) to 95°F (35°C)
Operating Relative Humidity Range	15% to 90%, non-condensing
Operating Atmospheric Pressure Range	1013 hPa to 715 hPa 0-3000 m (0 ft - 9840 ft)
Power Consumption	290 Watts Avg at nominal rated AC input, 25°C, 935 hPa
Weight	30 lbs. (13.6 Kilograms)
Sound Pressure Level at 3 and 5 LPM (ISO 80601-2-69)	<46dBA
Sound Power Level at 3 and 5 LPM (ISO 80601-2-69)	<54dBA
Alarm Sound Level	> = 55 dBA
Dimensions	18.4"H x 14.6"W x 9.9"D (46.6 x 37.0 x 25.0 cm)
Maximum Limited Pressure	Normal Condition: 9.25 PSIG (63.8 kPa), Single Fault Condition: 45 PSIG (310.3 kPa)
Operating System	Time Cycle / Pressure Swing
Low Oxygen Indicator	<85% Yellow - Low Oxygen <60% Red - Service Required
Transport and Storage Temperature Range	-13°F (-25°C) to 158°F (70°C)
Transport and Storage Relative Humidity Range	15% to 90%, non-condensing

CATALOG NUMBER	555DS
Equipment Class and Type	<input type="checkbox"/> Class II Equipment;  Type BF Applied Part IP21 ingress protection
Approval Body and Safety Standards	TUV IEC 60601-1 IEC 60601-1-6 IEC 60601-1-11 ISO 80601-2-69
EMC Compliance	EN 60601-1-2

Specifications subject to change without notice.

Units of measure for volumes, flows and leakages are expressed in STPD: 101.3 kPa at an operating temperature of 20 °C, dry Measurement Uncertainty

- Flow +/- 1.25% of reading or 0.0005 SLPM
- Pressure +/- 50 PPM of reading
- Oxygen purity +/- 0.1% oxygen

### OXYGEN CONCENTRATION VS FLOW RATE (Across the listed voltage)

555DS	
Flow L/m	%O <sub>2</sub>
5	90.6%
4	94.5%
3	94.8%
2	93.9%
1	93.6%
0.5	91.1%

**NOTE** - The %O<sub>2</sub> values were obtained through verification testing of the device's performance, and that the actual value may differ slightly from what is reported in the manual depending on the age of the device and environmental conditions of use.

### ELECTROMAGNETIC COMPATIBILITY INFORMATION

IMMUNITY TEST LEVELS			
TEST DESCRIPTION	SPECIFICATION	NOTES	RESULTS
CISPR11 & IEC 60601-1-2:2014 + AMD1:2020 (EDITION 4.1)			
RF Emissions	CISPR 11	Group 1	Complies
RF Emissions	CISPR 11	Class B	Complies
Harmonic Current Emissions	EN 61000-3-2:2014	AC Input Class B	Complies
Voltage Fluctuations & Flicker Emissions	EN 61000-3-3:2013	AC Input	Complies
Electrostatic Discharge Immunity (ESD)	EN 61000-4-2:2008	±15kV (Air) ±8kV (Contact)	Complies
Radiated Electromagnetic Field Immunity	EN 61000-4-3:2006	80MHz to 2.7GHz @ 10V/m 2.7-6.0GHz @ 3V/m 80% AM @ 1kHz	Complies
Fast Transient/Burst Immunity	EN 61000-4-4:2004	±2kV for Power Lines ±1kV for I/O Lines	Complies
Surge Immunity	EN 61000-4-5:2006	±1kV differential ±2kV common	Complies
Conducted RF Immunity	EN 61000-4-6:2009	3Vrms 150kHz to 80MHz 6Vrms ISM Band frequencies between 150kHz-80MHz	Complies

Maximum length of AC line cord = 3m (9.8 ft)

## **CYBERSECURITY STATEMENT**

The 555 Series 5L Compact Oxygen Concentrator is designed with limited connectivity for enhanced security. This device has no provision for network connection.

Do not use if the device screws or labels show signs of tampering. Contact Drive DeVilbiss at Consumersupport@drivemedical.com.

Please contact Drive DeVilbiss at Consumersupport@drivemedical.com for a copy of the Software Bill of Materials.

## **WARRANTY**

DeVilbiss Healthcare warrants the DeVilbiss 5 Liter Oxygen Concentrator under the conditions and limitations stated below. DeVilbiss warrants this equipment to be free from defects in workmanship and materials for three (3) years from date of factory shipment to the original purchaser, (typically the healthcare provider) unless contractually specified otherwise. This warranty is limited to the Buyer of new equipment purchased directly from Drive DeVilbiss, or one of its Providers, Distributors, or Agents. DeVilbiss' obligation under this warranty is limited to product repair (parts and labor) at its factory or at an Authorized Service Center. Routine maintenance items, such as filters, are not covered under this warranty, nor does it cover normal wear and tear.

## **WARRANTY CLAIMS SUBMISSION**

The original purchaser must submit any warranty claim to Drive DeVilbiss or to an Authorized Service Center. Upon verification of the warranty status, instructions will be issued. For all returns, the original purchaser must (1) properly package the unit in a DeVilbiss approved shipping container, (2) properly identify the claim with the Return Authorization Number, and (3) send the shipment freight prepaid. Service under this warranty must be performed by DeVilbiss and/or an Authorized Service Center.

**i NOTE** – This warranty does not obligate DeVilbiss to provide a loaner unit during the time that an oxygen concentrator is undergoing repair.

**i NOTE** – Replacement components are warranted for the unexpired portion of the original Limited Warranty.

- This warranty shall be voided, and DeVilbiss shall be relieved of any obligation or liability if:
- The device has been misused, abused, tampered with, or used improperly during this period
- Malfunction results from inadequate cleaning or failure to follow the instructions
- The equipment is operated or maintained outside the parameters indicated in the DeVilbiss operating and service instructions
- Unqualified service personnel conduct routine maintenance or servicing.
- Unauthorized parts or components (i.e., regenerated sieve material) are used to repair or alter the equipment.
- Unapproved filters are used with the unit.

THERE IS NO OTHER EXPRESS WARRANTY IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state.

**i NOTE** – International warranties may vary.

## **ORDERING AND RETURNING PARTS**

### **DeVilbiss Customer Service Contact Information**

Customer Service (USA): 800-338-1988

**International Department: 814-443-4881 / DHCinternational@DeVilbissHC.com**

Ordering Non-Warranty Replacement Parts

Order non-warranty parts and literature from your DeVilbiss provider

## **RETURN AND DISPOSAL**

This device may not be disposed of with household waste. After use of the device, please return the device to the provider for disposal. This device contains electrical and/or electronic components that must be recycled per EU Directive 2012/19 EU-Waste Electrical and Electronic Equipment (WEEE). Non-infectious used accessories (e.g. nasal cannula) can be disposed of as residential waste. The disposal of infectious accessories (e.g. nasal cannula from an infected user) must be made via an approved waste disposal company. Names and addresses can be obtained from the local municipality.

Store used sieve beds in sealed airtight container to minimize release of adsorbed substances. Waste must be handled in accordance with all federal, state, provincial, and local regulations. There are no known hazardous contaminants adsorbed into the molecular sieve therefore used sieve beds may be disposed of as residential waste. Never dispose by means of public sewers or drainage.

## **PROVIDER'S NOTES - Cleaning and Disinfection When There is a Patient Change**

DeVilbiss Healthcare recommends that at least the following procedures be carried out by the manufacturer or a qualified third party between uses by different patients.

The healthcare provider is to assess the needs of the patient, at time of installation, for backup supplies of supplementary oxygen in case of oxygen concentrator or power failure, based on the condition of the patient, the environment in which the patient lives and the ability to resupply the patient with backup supplies of supplementary oxygen. The healthcare provider should reassess as conditions change.

	<b>DANGER</b> Before attempting any cleaning procedures, turn the unit OFF and disconnect from mains power. Risk of electrical shock.
--	--

- i** **NOTE** – Recommendations for preventive maintenance at 3-year intervals are outlined in the Service and Maintenance Guidance below.
  - i** **NOTE** – If the following described complete processing of the concentrator by an appropriately trained individual is not possible, the device should not be used by another patient.
  - i** **NOTE** – If preventive maintenance is due at this time, these procedures should be carried out in addition to the servicing procedures.
1. Use disinfectants safely. Always read the label and product information before use.
  2. Always wear personal protective equipment when performing this procedure. Use suitable gloves and safety glasses. Cover exposed skin on arms to prevent accidental contact with bleach solution that has been applied to the concentrator.
  3. Dispose of all accessories that are not suitable for reuse. This includes but may not be limited to the oxygen tubing, tubing connectors, nasal cannula and/or mask, oxygen outlet connector, and humidifier bottle. Reusing accessories marked for single use may cause risk of cross contamination or loss of performance.
  4. Clean the exterior of the concentrator with a clean lint-free cloth. Heavy soil should be removed with a clean lint-free cloth dampened with clean tap water. A soft bristled brush dampened with clean tap water can be used to remove stubborn soil. Dry the concentrator using a clean lint-free cloth if clean tap water was used to remove soil.
  5. Use 5.25% chlorine bleach (Clorox Regular Liquid Bleach or equivalent). Mix one (1) part bleach with four (4) parts clean tap water in an appropriate clean container. This ratio produces a one (1) part bleach to five (5) total parts solution (1:5). The total volume (amount) of solution required is determined by the number of concentrators in need of disinfection.
  6. Apply the bleach solution in an even manner to the exterior cabinet and power cord using a clean lint-free cloth. The cloth should be dampened only and not dripping of solution. Do not use a spray bottle to apply the solution. Do not saturate the device with the solution. Take care that no solution enters the vent areas on the concentrator base. Avoid over-saturating the exterior cabinet seams so that no solution residue builds up in these areas. Avoid the caster wells located on the bottom of the unit.
  7. Exposure time of the disinfectant solution should be 10 minutes minimum to 15 minutes maximum.
  8. After the recommended exposure time, all surfaces of the concentrator should be wiped with a clean lint-free cloth dampened with drinking quality water no warmer than room temperature. Dry the unit with a dry, clean lint-free cloth. This is to remove residue that may stain or leave a film on the unit, especially after repeated disinfections.
  9. Check the power cord, the plug on the back of the device, the power switch, the circuit breaker, and the LCD display for possible damage. Replace all damaged or worn components.
  10. Check the oxygen concentration. If the device is within specification, the extended life intake bacteria filter does not need to be replaced between patients. If the oxygen concentration is not within specification, the provider should refer to the service manual section on Troubleshooting.

- i NOTE** – There is no portion of the gas pathways through the concentrator that could be contaminated with body fluids under normal conditions.
- The device patient connection may unintentionally become contaminated with expired gases for a single fault condition i.e., a hose internal to the device becomes disconnected. This condition will cause no flow out of the device and/or an alarm condition. Should this occur, refer to the service manual for additional instructions.

## **DISINFECTION**

- i NOTE** – The disinfection process can only be completed by the manufacturer or by an appropriately trained service technician.

	<b>Recommended disinfection interval</b>	<b>Number of disinfection cycles</b>	<b>Compatible disinfection method</b>
Exterior Cabinet, power cord	Between patients	20	1:5 chlorine bleach (5.25%) and water solution
Oxygen tubing, tubing connectors, nasal cannula/mask, oxygen outlet connector, humidifier bottle	Do not clean, replace between patients	N/A	N/A

- i NOTE** – Frequency of inspection and cleaning of filter may be dependent upon environmental conditions like dust and lint. The recommended disinfection interval may be performed when a device is used on the same patient, as determined by clinical staff. Disposable accessories should be replaced according to the accessory manufacturer's instructions.

## **SERVICE AND MAINTENANCE GUIDANCE**

Service and maintenance should only be performed by appropriately trained and authorized Drive DeVilbiss personnel and/or service centers.

<b>DeVilbiss Oxygen Concentrator Preventive Maintenance/Service Guide</b>								
<b>Model</b>	<b>Oxygen Purity Verification</b>	<b>Intake HEPA Filter</b>	<b>Internal Compressor Filter</b>	<b>Final HEPA Filter *</b>	<b>Sieve Beds *</b>	<b>Compressor Cup Seals *</b>	<b>Audible and visual Indicators</b>	<b>AC power cord, power switch, and circuit breaker</b>
555 Series	Every 3-years or between patient uses, whichever comes first.	Inspect every 3- years or between patient uses, whichever comes first. Replace if filter is dirty or damaged or there is a drop in O2 purity.	Inspect in conjunction with compressor service. Replace if filter is dirty or damaged.	Inspect in conjunction with compressor service. Replace if filter is dirty or damaged	When indicated by device performance below specification for oxygen purity, operating pressures and/ or other indications of component wear.	When indicated by device performance below specification for oxygen purity, operating pressures and/or other indications of component wear.	Inspect every 3 years. Replace if defective or damaged.	Inspect between patients. Replace if defective or damaged.

\* Sieve bed, compressor cup seal, compressor filter and final HEPA filter service should only be performed by appropriately trained and certified Drive DeVilbiss service centers.

- i NOTE** – This is a suggested maintenance and service schedule for home oxygen providers. Individual maintenance requirements may vary based upon local operating conditions, regulations, or other circumstances.

## **INITIAL INSPECTION**

- Upon receiving, examine the unit for external damage. If the unit appears to have external damage, please contact DeVilbiss for assistance.
- Check to be sure the intake filter is in place.
- Plug the unit into an electrical outlet, turn the unit "ON" and check the audible/visual alarms. When the unit is turned "ON," the LCD display shows a Drive DeVilbiss logo and the hour meter and an audible signal will briefly alarm confirming that the LCD display and audible signal are functioning properly. The unit will then operate in "start up" mode with the Startup Period Low O2 symbol shown on the LCD display until a normal oxygen level is achieved. When a normal oxygen level is achieved, the Startup Period Low O2 symbol turns off and the Green Normal Oxygen symbol shows on the LCD display. The "start up" period may take up to 15 minutes. If the LCD display DOES NOT illuminate or the audible alarm DOES NOT sound, this indicates the alert system is not functioning properly. Refer to the Troubleshooting chart or contact DeVilbiss for assistance.
- Set the oxygen flow at the maximum recommended flow rate and allow the unit to run for 20 minutes. The internal oxygen sensor monitors the oxygen purity. If the oxygen is within specification, the Green Normal Oxygen symbol will be illuminated. If the Yellow Low O2 symbol is illuminated, refer to the Service Manual or contact DeVilbiss for assistance.

5. With unit still running, unplug to test the power fail alarm. If the power fail alarm does not provide an audible alert, refer to the Service Manual or contact DeVilbiss for assistance.

**i** **NOTE:** *When there is a power outage, turn the unit's power switch to OFF. When power is restored, turn the unit's power switch to ON to check the audible/visual alarm functions.*

## **OXYGEN PROVIDER PREVENTIVE MAINTENANCE GUIDANCE**

**i** **NOTE** – *Scheduled maintenance should be performed in accordance with the Preventive Maintenance/Service Guide table above or Between Patient Uses.*

1. Discard all oxygen tubing, cannula/mask, oxygen outlet connector and humidifier bottle.
2. Follow the Cleaning and Disinfection Instructions in the IFU.
3. Clean the concentrator exterior cabinet and inspect/replace filters in accordance with the table above.
4. Inspect all plugs, cords, and components. Replace any damaged or worn components.
5. Check oxygen concentration with a calibrated oxygen analyzer and record the oxygen percentage. If the concentration is not within specification, refer to troubleshooting section of the IFU or the Service Manual.
6. Record the unit hours of use, the hour meter is shown on the LCD display for a few seconds upon initial power up.
7. Verify Audible Alert and LCD display at each service at startup and while operating.
8. Turn the oxygen concentrator on and allow the unit to run for at least 20 minutes. With the unit still running, unplug to test the power fail alarm. If the power fail alarm does not provide an audible alert, refer to the Service Manual or contact DeVilbiss for assistance.

Install, use and maintain the 555 Series 5L Compact Oxygen Concentrator following this instruction guide to minimize the environmental impact of the oxygen concentrator through its expected life. The oxygen concentrator consumes 290 watts maximum during normal use. Distilled water can be used as needed for humidification of the delivered oxygen. The Intake Filter and Final Bacterial Filter are disposable and could need to be replaced twice during the expected service life. Other byproducts created during the normal use of the oxygen concentrator are heat, noise and nitrogen gas. Follow instructions to minimize the effects of heat and noise. Nitrogen gas discharged by the oxygen concentrator is quickly dissipated back into room air. The oxygen concentrator contains no hazardous substances.

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
**Drive DeVilbiss Healthcare**

Corporate Headquarters  
99 Seaview Blvd  
Port Washington, NY 11050  
USA  
877-224-0946 • 516-998-4600  
[www.drivemedical.com](http://www.drivemedical.com)

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A-555DS Rev D

©2026 Medical Depot, Inc. dba Drive Medical, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive Medical are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies.

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**DeVilbiss Healthcare LLC**