



## Instructions for Placing the iGo®2 Portable Oxygen Concentrator (POC) in the Backpack

### LOWER COMPARTMENT:

1. Locate the two padded legs inside the lower compartment of the backpack (Figure 1). Extend each leg downward (Figure 2). Attach and secure each leg to the Velcro strips on the left and right sides of the lower compartment (Figure 3). Securing these legs will provide additional support for the lower compartment, maximizing the storage space.
2. Use the lower compartment of the backpack to store the AC power cable, DC power cable, and any extra batteries (if applicable).



Figure 1.



Figure 2.



Figure 3.

### UPPER COMPARTMENT:

1. Place the iGo®2 Portable Oxygen Concentrator (POC) into the upper compartment of the backpack.
2. Align with the oxygen outlet port (Figure 4) and power input (Figure 5), then secure the iGo®2 unit into the upper compartment using the two Velcro straps.



Figure 4.



Figure 5.

**WARNING:** Make sure the air inlet and exhaust vents are not blocked. The exhaust vents on the iGo®2 unit will align at the bottom of the backpack. These vents can emit hot exhaust gas. Do not place near bare skin.

### WEARING THE IGO®2 BACKPACK

1. Slide your arms through the straps and align the iGo®2 backpack on your back. Adjust the straps to a position where the iGo®2 backpack is comfortable.