



**ResMed**

**CPAPXCHANGE**

SLEEP WELL TONIGHT



**AirTouch™** F20

Full face mask

[ResMed.com](https://www.ResMed.com)

# Start dreaming again

The benefits of sleep apnea treatment are truly remarkable. Getting a good night's sleep gives you the energy to get back to doing the things you love. But one of the most inspiring things people tell us about starting CPAP is that they're able to dream again.






Congratulations on taking an important step toward getting restful sleep.

ResMed is the #1 mask brand preferred by patients,<sup>1</sup> and we're excited to present the AirTouch™ F20 – the softest mask we've ever designed.

Featuring a ResMed UltraSoft™ memory foam cushion that provides a uniquely comfortable and personalized fit, the AirTouch F20 is a mask we're confident will help you get the good night's sleep you deserve.

**AirTouch™** F20  
Full face mask

## ResMed UltraSoft memory foam for superior sleep



The AirTouch F20 cushion is made of UltraSoft memory foam for a light, breathable fit and seal. A simple wipe down keeps the cushion clean, and regular replacement keeps it maintenance-free.

## Changing the face of comfort

The AirTouch F20's UltraSoft memory foam cushion adapts to the unique contours of your face, creating a seal that is designed to work across a range of therapy pressures. The AirTouch F20 has also proven to work for nearly everyone who tried it – fitting 98% of patients in a recent study.<sup>2</sup>



### **Magnetic clips**

guide the headgear to the frame quickly and easily.

### **Quick-release elbow**

simple to disconnect from the mask tubing without removing the mask.



# The importance of part replacement

Learn why regular mask part replacement is especially important for the AirTouch F20's UltraSoft memory foam cushion.

Register for replacement reminders at [ResMed.com/AirTouchReminders](https://www.ResMed.com/AirTouchReminders)



**1** Combined responses of those expressing a preference in an independent patient survey in the UK, USA, France, Germany and Australia between April 10 and May 11, 2015. Visit [ResMed.com/MaskBrand](https://www.ResMed.com/MaskBrand). **2** ResMed AirTouch F20 internal USA fitting study of new and existing CPAP patients, conducted April 2016 and August 2016.

Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA  
+1 858 836 5000 or 1 800 424 0737 (toll free). See [ResMed.com](https://www.ResMed.com) for other ResMed locations worldwide. AirTouch and UltraSoft are trademarks and/or registered trademarks of the ResMed family of companies. Specifications may change without notice. © 2017 ResMed. 1019465/1 2017-01

[ResMed.com](https://www.ResMed.com)

Follow ResMed on:

