


Printing instructions

1 Print the sizing guide

- Print on U.S. Letter- or A4-sized paper.
- Ensure you select "Actual size" and not "Fit" in the printer settings.
- For best results, print in color.




2 Measure your printout

- Place your ruler on the scale below and ensure they measure the same.

Sizing guide instructions


1 Carefully cut along the dotted lines.




2 Hold the guide so that it presses against the underside of your nose.

TIP

Use a mirror or ask someone at home to assist you with sizing.



CPAPXCHANGE
SLEEP WELL TONIGHT



3 Choose the size that best fits your nose.