

F&P Nova Nasal



Nasal Mask

SwingFit™.
Performance defined.



Fisher & Paykel
HEALTHCARE

Discover the definition of performance

For a sound sleep, the Mask Matters Most™. Introducing F&P Nova™ Nasal with SwingFit headgear. From optimal fitting and seal to minimal noise and draft – it's performance defined.



Performance you can trust

98%

of trial users trust F&P Nova Nasal to work¹

Easy fitting

Sound sleep

Stable seal

Combining innovative technologies


Packed with features working in synergy with the intuitive SwingFit headgear, F&P Nova Nasal delivers all-around performance.

SwingFit headgear

Simple one-click fitting with a clear field of view.



Tube clip
Optional over-the-head tube placement.



RollFit™ cushion with pivoting clips

Adaptive movement for a stable seal in all directions.

Diffuser

Minimizes noise and draft.

Scan to find out more.




F&P Nova Nasal

Product Codes

Code	Description
NVNISA	Nova Nasal Mask - Small
NVNIMA	Nova Nasal Mask - Medium
NVNILA	Nova Nasal Mask - Large
NVNISMLA	Nova Nasal Mask - Fit pack - Small, Medium, Large

Spare Parts

Code	Description
400NVN114	Nova Nasal Spare Cushion - Small
400NVN115	Nova Nasal Spare Cushion - Medium
400NVN116	Nova Nasal Spare Cushion - Large
400NVN125	Nova Nasal Headgear - Standard
400NVN126	Nova Nasal Headgear - Large
400NVN131	Nova Nasal Mask No Headgear - Small
400NVN132	Nova Nasal Mask No Headgear - Medium
400NVN133	Nova Nasal Mask No Headgear - Large
400NVN123	Nova Nasal Back Panel - Standard*
400NVN124	Nova Nasal Back Panel - Large*
400NVN141	Nova Nasal Frame Assembly Spare
400NVN151	Nova Nasal Diffuser Spare
400NVN161	Nova Nasal Tube Clip Spare
400NVN162	Nova Nasal Headgear Pad Spare




F&P SureFit
30-Day
 Mask Guarantee



F&P myMask™ app

Supports your patients with mask selection, sizing and step-by-step videos to assist their mask setup and use. Scan the QR code or download directly.



* Headgear includes back panel. Back panel can also be purchased separately.

Always read the label and follow the directions for use. Prescription only. This product may not be right for you. Your healthcare provider will advise you whether this product is suitable for you/your condition. If symptoms persist, talk to your healthcare provider.

1.42 out of 43 participants agreed or strongly agreed with the following statement: "I trust the F&P Nova Nasal mask to work." Internal validation trial conducted by Fisher & Paykel Healthcare on 43 participants in the United States, 2023.

F&P, F&P Nova, Nova, Mask Matters Most, myMask, SureFit, SwingFit and RollFit and are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. The App Store is a service mark of Apple Inc., in the U.S. and other countries. Google and Google Play are trademarks of Google LLC.

633609 REV A © 2025 Fisher & Paykel Healthcare Limited.

Fisher & Paykel
 HEALTHCARE