

The logo consists of the text "CPAPX" in a bold, sans-serif font. The "CPAP" is in white, and the "X" is in blue. The text is set against a dark blue background that is shaped like a right-pointing arrow.

Transforming sleep.

Forma™

FULL FACE MASK

The Forma™ Full Face Mask forms a superior seal and provides even more comfort with a greater range of movement through the active contouring of its new FlexiFoam™ Cushion. The T-Piece adds stability to the stylized Mask Base, and works in harmony with the highly contoured Silicone Seal and Under-Chin Design to enhance sleep performance. With its ease of use and superior conforming seal, the Forma™ sets a whole new standard.


Transform your patient's sleep experience today, with the Fisher & Paykel Healthcare Forma™ Full Face Mask.

FlexiFoam™ Cushion & T-Piece Mask Base




Fisher & Paykel
HEALTHCARE

Full Face Delivery – *Forma™*




The Forma™ Full Face Mask offers an improved seal and a more comfortable fit due to the interaction between the new FlexiFoam™ Cushion and the Under-Chin Design. This alternative to nasal or oral delivery comes fully assembled and ready for use.

UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




FlexiFit™ Technology Auto-contouring fit

The FlexiFoam™ Cushion using FlexiFit™ Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments.




FlexiFoam™ Cushion Active-contouring cushion

The FlexiFoam™ Cushion provides greater comfort with the unique design of the soft, light and supportive foam. The FlexiFoam™ features the advanced active contouring that conforms naturally to your patient's face.




Glider™ Strap with Easy-Release Cord

The sliding Glider™ Strap with Easy-Release Cord provides freedom of movement and one-step removal. The mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks. Simply pull the Easy-Release Cord for one-step removal.




Under-Chin Design with Crown Strap

An Under-Chin fit provides extra stability and maintains the seal, even with movement of the jaw during the night. The Crown Strap works in conjunction with the Under-Chin Design to eliminate the need for a chin strap.




T-Piece Mask Base Optimal seal

With no complicated adjustments, the T-Piece provides added stability and complements the Under-Chin Design to allow FlexiFit™ to auto-adjust for hassle-free fitting.



Advanced Air Diffuser Quiet usage

Maintenance Free. This efficient Air Diffuser system prevents air from blowing onto the sleeping partner and decreases noise that could potentially disrupt sleep.



Frosted Silicone Seal Optimizing seal

The enhanced seal technology has a soft, smooth feel and slides easily onto the face. Its frosted texture optimizes mask seal performance.

CONTENTS

- Complete mask with headgear assembled and ready for use
- Four sizes: Small, Medium and Large (two sizes included) and X-Large, sold separately
- User-friendly Fitting Guide

REORDER CODE:

400470 Small
400471 Medium and Large
400473 X-Large

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

Manufacturer Fisher & Paykel Healthcare Limited, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand
Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

Australia Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Northern Ireland** Tel: 0800 132 189 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Fisher & Paykel Healthcare UK Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA** Tel: 1 800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001