


Handheld Mesh Nebulizer with Rechargeable Battery

Item#
NEB402



- Includes internal rechargeable battery that provides up to 1 hour of continuous use after full charge
- Great value pricing and quality
- Easy to operate and transport

Included

- Medication cup
- USB cord
- Mask
- Mouthpiece
- Carry bag
- Manual



Small particle size

Nebulizes a fine mist of <5µm particles to further reach and treat respiratory system.

Quick treatment

Saves time with a faster nebulization rate than standard nebulizers.

Lightweight and discreet

Weighs just 1/4 lb for greater portability and features quiet operation.

Particle Size	MMAD < 5µm
Nebulization Rate	0.25ml/min-0.9ml/min
Medication Cup Capacity	10ml (Max)
Device Weight	Approx. 1/4 lb
Device Size	Approx. 1 1/2 x 2 x 3 7/8 inches