

Nasal Aire II[®] Headgear

Multiple Wearing Options

Standard Method :

All components included in Package



Use headband , two (2) tubing holder rings , neckstrap and bolo .

Florida Method :

* Optional headgear



Use Headband , two (2) figure 8 loops and bolo.

Ponytail Method :

All components included in Package



Use Headband , two (2) tubing holder rings and bolo.

Breatheband Method :

* Optional headgear



Use Breatheband and neckstrap.

Chinstrap Method :

* Optional headgear



Use chinstrap and neckstrap

Figure 8 Method :

All components packaged



Use figure 8 loop

* Denotes optional item. Sold Separately.