

**Proven  
best-in-class  
reliability,  
year after year.**



# Oxlife INDEPENDENCE<sup>®</sup>

## Portable Oxygen Concentrator

Our flagship device, the Oxlife **INDEPENDENCE**, sets the standard for **long operating life, reliability, durability and 24/7 performance.**

Equipped with DNA Technology, it is the **DME choice for hospital discharge** and the hospital-to-home transition, meeting the majority of long-term oxygen therapy requirements.



**Continuous flow up to 3 LPM**  
and 12 pulse settings



**Proven 24/7 use, year after year**  
backed by a five-year warranty




**Use as POC & stationary**  
dual billing (E1390 & E1392)



**No tanks to refill**  
or ongoing tank deliveries



**24/7 remote monitoring**  
with DNA Technology




**Long-life, dual battery**  
operation of up to six hours



**Intuitive, easy to use**  
LCD display and setting controls





**DNA TECHNOLOGY®**

# Remote asset & clinical monitoring

O<sub>2</sub> Concepts' proprietary, device-initiated connectivity solution

---

GPS-enabled inventory control  
and loss prevention technology

---

24/7 remote monitoring  
for device performance

---

Easily access data from any device  
via secure, cloud-based platform

---

Proactively promote enhanced  
patient therapy adherence

---

Reduce field calls with intuitive  
customer support resources

---

Manage hospital DME closet  
inventory with device Geofence



## OXLIFE INDEPENDENCE® PRODUCT SPECIFICATIONS<sup>1</sup>

<b>Dimensions</b>	Height: 20.2" • Width: 10.9" • Depth: 9.45"
<b>Weight<sup>2</sup></b>	Device: 16.7 lbs
<b>Oxygen Flow Rate</b> (Pulse 15-40 BPM)	Continuous: 0.5 - 3 LPM in 0.5 increments
	Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96 ML)
<b>Battery Duration<sup>3</sup></b>	5.75 hours @ 2 pulse at 20 bpm
	2.5 hours @ 2 LPM continuous
<b>Battery Recharge Time<sup>4</sup></b>	2.5 hours
<b>Oxygen Concentration</b>	87% - 95% at sea level at all settings
<b>AC Power</b>	Universal Input: 100 - 240V, 50/60Hz
	50 - 60 Hz • 2.0A maximum
<b>DC Power</b>	12 - 15 VDC • 15 amps

<sup>1</sup> O<sub>2</sub> Concepts® reserves the right to revise specifications at any time without notice <sup>2</sup> Production weight excludes cart kit

<sup>3</sup> Battery duration presented for two-battery system <sup>4</sup> Battery charge time increases if device is operating while charging

## Smarter, more efficient and more reliable oxygen therapy than ever before

A legacy of design innovations, cutting edge technology and unparalleled support

Relied upon by leading academic and research medical centers, and by hospital groups and homecare providers nationwide

The most cost-effective device & operating model

A simplified, effective hospital/clinic discharge process

The ability to remotely monitor and locate every device in your fleet

Designed for the reliability and durability needs of the DME industry

Contact us at **(877) 867-4008 (prompt 6)** or at **info@o2-concepts.com** for more information or to consult with one of our territory managers.

**o2-concepts.com**



O<sub>2</sub> Concepts®, Oxlife INDEPENDENCE®, DNA Technology® and Energy Smart Technology® are trademarks of O<sub>2</sub> Concepts®, LLC and may not be reproduced without express consent.  
© 2024 O<sub>2</sub> Concepts®, LLC. All rights reserved.  
900-1019 Rev A