

**CPAPX**



**GOING GREEN TODAY**

**FOR A BETTER TOMORROW**



**AEROSOL DELIVERY SYSTEM**



PARI Vios® delivers consistent particle size with targeted aerosol delivery to the lungs in just 5-6 minutes.<sup>1</sup> Available in both standard and pediatric configurations, PARI Vios includes environmentally friendly packaging, instructions, and 100% recycled cotton tote to help reduce our carbon footprint.

[www.pari.com/vios](http://www.pari.com/vios)

# ADULT AND PEDIATRIC OPTIONS



Standard Configuration  
(Includes Back Up Nebulizer: not shown)



Pediatric configuration also includes  
Bubbles the Fish™ II mask, carry tote, and activity set

## Configurations

	Part #	HCPCS #
Standard, LC® Sprint	310F35-LCS	E0570, A7005
Standard, LC PLUS®	310F83-LC+	E0570, A7005
Pediatric, LC® Sprint	310F35-P	E0570, A7005, A7015
Pediatric, LC PLUS®	310F83-P	E0570, A7005, A7015

## Technical Specifications

Operating Pressure <sup>2</sup>	18 PSI (1.2 bar)
Liter Flow <sup>2</sup>	4.5 LPM
Weight	3 lbs
Dimensions (L, W, H)	6.5 x 6.5 x 3.75 in
Compressor Warranty	5 years

<sup>1</sup> Measured with Malvern MasterSizer X at 50% relative humidity, 0.9% NaCl solution, inspiratory flow 20 liters/minute, continuous nebulization, 23° C, fill volume 2.5 ml.

<sup>2</sup> When used with a PARI reusable nebulizer.

## PATIENTS BENEFIT FROM:



Fast Treatments

Fast 5-6 minute treatments help improve patients' quality of life



Breath-enhanced nebulizer design saves money by wasting less medication



Consistent particle size with targeted medication delivery to the lungs

## Accessories



PARI LC Sprint  
023F35



Bubbles the Fish™ II  
Pediatric Mask  
044F7248



Replacement  
Air Filter  
041F4851P2



PARI LC PLUS  
022F81



PARI Adult Aerosol Mask  
044F7252



**PARI** Respiratory Equipment, Inc.  
Tel: 1-800-FAST-NEB (327-8632)  
product.info@pari.com  
www.pari.com | Midlothian, VA 23112