

CPAPX



PARI VORTEX®

Antistatic Holding Chamber/Spacer

WITH NEW PEDIATRIC MASKS



www.pari.com/vortex

Advantages of PARI VORTEX®



- Antistatic, aluminum chamber ensures efficient and consistent medication delivery with MDIs for optimal treatments.
- The flexible connecting ring is compatible with most common metered-dose inhalers.
- Mouthpiece cap protects the VORTEX® from unwanted debris.
- The VORTEX® and masks are dishwasher safe and boilable, making cleaning and disinfection super easy!
- VORTEX grows with you—our interchangeable masks eliminate the need to purchase a new spacer when changes in mask size are needed.
- The patented one-handed operation aid makes treatments easier with babies or for those who have manual coordination difficulties.

PARI VORTEX® – ideal for patients of all ages



VORTEX with Pediatric Mask, Sm
0 - 1.5 years
Part No.: 051F8100
HCPCS No.: A4627



VORTEX with Pediatric Mask, Med
1 - 4 years
Part No.: 051F8200
HCPCS No.: A4627



VORTEX with Adult Mask
4 years and up
Part No.: 051F8300
HCPCS No.: A4627



VORTEX Pediatric Mask, Sm
0 - 1.5 years
Part No.: 044F8100
HCPCS No.: A7015



VORTEX Pediatric Mask, Med
1 - 4 years
Part No.: 044F8200
HCPCS No.: A7015



VORTEX with Mouthpiece
4 years and up
Part No.: 051F8000
HCPCS No.: A4627




One-Handed Operation Aid
Part No.: 044F5100

EXCELLENT AEROSOL PERFORMANCE

A good seal is key to effective inhalation when using a mask with a holding chamber/spacer. PARI's new masks show excellent sealing capabilities even when applied gently to the face. When compared to MDI alone, PARI Vortex with masks show an increased fine particle dose (FPD, $<5\mu\text{m}$) and low throat deposition with both salbutamol and budesonide.

Figure 1: Results for FPD of pMDI, VORTEX® with baby mask and VORTEX® with child mask for salbutamol and budesonide



DISHWASHER SAFE


Thanks to the durable aluminum chamber and removable mouthpiece, the VORTEX® and masks can be boiled, disinfected, and washed in the dishwasher.

PARI LENDS A HAND WITH THE ONE-HANDED OPERATION AID

The patented one-handed operation aid helps simplify treatments with babies and small children as well as for adults who may have manual coordination issues or reduced strength in their fingers.



VORTEX® pediatric masks



Aluminum Chamber. Less Plastic. Better Delivery.



PARI Respiratory Equipment, Inc.
Tel: 1-800-FAST-NEB (327-8632)
product.info@pari.com
www.pari.com | Midlothian, VA 23112