

P2-E7 PORTABLE OXYGEN CONCENTRATOR



ITEM #	DESCRIPTION
P2-E7	7 Setting Portable Oxygen Concentrator

- Compact & Lightweight.
- High efficiency, low noise. Excellent Heat Dissipation.
- 93% O₂ Concentration.
- High Adsorption Sieve Bed.
- Includes AC Power Supply & DC Car Charger.
- Conforms to all applicable FAA compliant criteria for POC carriage and use on board aircraft.

FEATURES

Size	8.7" x 3.4" x 6.3"
Weight	4.37lbs
Warm up time	2 Minutes
Settings	Seven Flow Settings (1-7)
Maximum Oxygen Output	1400 ml/min
Oxygen Concentrator	87%-96%
AC Power Supply	100-240v, 50-60hz (auto sensing to allow for worldwide use)
DC Power Supply	Input: 11-16VDC, Output: 19V 6.3a
Rechargeable Battery	14.4VDC
Battery Rated Capacity	6.7Ah
Display Screen	2.8"
Battery	6700mAh, Duration up to 7 Hours on Setting 1. Approx. 4 Hours Charge Time with AC.
Operation	Simple Control Functions and Easy-to-Read LCD Display
Use	Designed for 24/7 use

Setting	Battery Life
1	7 h
2	5.5 h
3	3.7 h
4	3 h
5	2.5 h
6	2 h
7	1.7 h

WARRANTY

Device: 5 Years
 Sieve Bed: 2 Years
 Battery: 1 Year

Applies to units purchased on or after January 1, 2026. Prior purchases remain subject to original warranty terms.

RHC-POC-005: V12: 2/4/26