

Rio™ II

Full Face Mask



The React Health Rio™ II Full Face Mask is a minimal contact mask that minimizes nose discomfort and doesn't interfere with vision or glasses.



Rio™ II
Full Face Mask

- ✦ Under-The-Nose Cushion for enhanced comfort by freeing the nasal bridge, which may reduce feelings of claustrophobia.
- ✦ Circular dispersion ventilation system evenly distributes airflow while minimizing noise- as low as 17 Db(A)
- ✦ Modular main frame compatible with all cushion sizes. (S,M,L)
- ✦ Adjustable headgear that provides a comfortable, secure fit for a wide range of head sizes.
- ✦ Features easy release clips- not magnets - that are common to all React Health masks.

Rio™ II Full Face Mask



Equipped with innovative magnet-free clips



Circular dispersion ventilation system



Designed for effortless disassembly, facilitating thorough and efficient cleaning



www.reacthealth.com info@reacthealth.com
Phone: 863-226-6285 Fax: 863-226-6284



MASK PART NUMBERS

RFF1000 - Rio II Full Face Mask - All Size Kit

RFF1001 - Rio II Full Face Mask - SM

RFF1002 - Rio II Full Face Mask -MD

RFF1003 - Rio II Full Face Mask - LG

REPLACEMENT SUPPLIES PART NUMBERS

RFF2001 - Rio II Full Face Repl Headgear

RFF3001 - Rio II Full Face Repl Cushion - SM

RFF3002 - Rio II Full Face Repl Cushion - MD

RFF3003 -Rio II Full Face Repl Cushion- LG

RFF4001 - Rio II Full Face Mask w/o Headgear - SM

RFF4002 - Rio II Full Face Mask w/o Headgear - MD

RFF4003 - Rio II Full Face Mask w/o Headgear - LG

ReactHealth.com

© 2024, React Health. All rights reserved.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change.

MI.MK.0392, Rev A