



SILENT NIGHT FULL FACE MASK COMPATIBILITY (UPDATED 2025-0109)

BRAND	MASK	SIZE	BEST LINER OPTION	OVERALL FIT	NOTES
Fisher & Paykel	Flexifit 432	S	Large (F103)	Great	The lower Quick Stick tabs of the Large liner will be able to reach the hard plastic on the bottom side of the mask. Even a small amount of adhesive contact to the hard plastic is adequate to hold the liner in place.
Fisher & Paykel	Flexifit 432	M	Large (F103)	Great	The lower Quick Stick tabs of the Large liner will be able to reach the hard plastic on the bottom side of the mask. Even a small amount of adhesive contact to the hard plastic is adequate to hold the liner in place.
Fisher & Paykel	Flexifit 432	L	Large (F103)	Great	The lower Quick Stick tabs of the Large liner will be able to reach the hard plastic on the bottom side of the mask. Even a small amount of adhesive contact to the hard plastic is adequate to hold the liner in place.
Fisher & Paykel	Forma	S	Large (F103)	Great	The lower Quick Stick tabs of the Large liner will be able to reach the hard plastic on the bottom side of the mask. Even a small amount of adhesive contact to the hard plastic is adequate to hold the liner in place.
Fisher & Paykel	Forma	M	Large (F103)	Great	The lower Quick Stick tabs of the Large liner will be able to reach the hard plastic on the bottom side of the mask. Even a small amount of adhesive contact to the hard plastic is adequate to hold the liner in place.
Fisher & Paykel	Forma	L	Large (F103)	Great	The lower Quick Stick tabs of the Large liner will be able to reach the hard plastic on the bottom side of the mask. Even a small amount of adhesive contact to the hard plastic is adequate to hold the liner in place.
Fisher & Paykel	Simplus	S	Small (F101)	Excellent	
Fisher & Paykel	Simplus	M	Medium (F102)	Excellent	
Fisher & Paykel	Simplus	L	Large (F103)	Excellent	
ResMed	AirFit F10	XS	Small (F101)	Excellent	
ResMed	AirFit F10	S	Small (F101)	Excellent	
ResMed	AirFit F10	M	Medium (F102)	Great	
ResMed	AirFit F10	L	Large (F103)	Great	
ResMed	AirFit F20	S	Small (F101)	Excellent	
ResMed	AirFit F20	M	Medium (F102)	Excellent	
ResMed	AirFit F20	L	Large (F103)	Excellent	
ResMed	AirFit F30	S	One-Size-Fits-All (V200)	Good	
ResMed	AirFit F30	M	One-Size-Fits-All (V200)	Good	
ResMed	AirTouch F20	S	Small (F101)	Excellent	Align the top edge of the liner with the top of the memory foam.
ResMed	AirTouch F20	M	Medium (F102)	Excellent	Align the top edge of the liner with the top of the memory foam.
ResMed	AirTouch F20	L	Large (F103)	Excellent	Align the top edge of the liner with the top of the memory foam.
ResMed	Mirage Quattro	XS	Small (F101)	Excellent	
ResMed	Mirage Quattro	S	Small (F101)	Excellent	
ResMed	Mirage Quattro	M	Large (F103)	Excellent	
ResMed	Mirage Quattro	L	Large (F103)	Excellent	
ResMed	Quattro Air	XS	Small (F101)	Excellent	Click here to see an FAQ with attachment options for the Quattro Air.
ResMed	Quattro Air	S	Small (F101)	Excellent	Click here to see an FAQ with attachment options for the Quattro Air.
ResMed	Quattro Air	M	Medium (F102)	Excellent	Click here to see an FAQ with attachment options for the Quattro Air.
ResMed	Quattro Air	L	Large (F103)	Excellent	Click here to see an FAQ with attachment options for the Quattro Air.
ResMed	Quattro FX	S	Small (F101)	Excellent	Attach upper Quick Stick tabs under head gear wings.
ResMed	Quattro FX	M	Medium (F102)	Excellent	Attach upper Quick Stick tabs under head gear wings.
ResMed	Quattro FX	L	Large (F103)	Excellent	Attach upper Quick Stick tabs under head gear wings.
ResMed	Ultra Mirage	S	Small (F101)	Excellent	
ResMed	Ultra Mirage	M	Medium (F102)	Excellent	
ResMed	Ultra Mirage	L	Large (F103)	Excellent	
Respironics	Amara Full Face	P	Small (F101)	Excellent	
Respironics	Amara Full Face	S	Small (F101)	Excellent	
Respironics	Amara Full Face	M	Medium (F102)	Excellent	
Respironics	Amara Full Face	L	Large (F103)	Excellent	
Respironics	Amara View	S	One-Size-Fits-All (V200)	Excellent	
Respironics	Amara View	M	One-Size-Fits-All (V200)	Excellent	
Respironics	Amara View	L	One-Size-Fits-All (V200)	Excellent	
Respironics	Comfort Gel Full	S	Small (F101)	Excellent	
Respironics	Comfort Gel Full	M	Medium (F102)	Excellent	
Respironics	Comfort Gel Full	L	Large (F103)	Excellent	
Respironics	DreamWear Full Face	S	One-Size-Fits-All (V200)	Excellent	Quick Stick tabs may only partially engage with the hard plastic.
Respironics	DreamWear Full Face	M	One-Size-Fits-All (V200)	Excellent	Quick Stick tabs may only partially engage with the hard plastic.
Respironics	DreamWear Full Face	MW	One-Size-Fits-All (V200)	Great	Quick Stick tabs may only partially engage with the hard plastic.
Respironics	DreamWear Full Face	L	One-Size-Fits-All (V200)	Great	Quick Stick tabs may only partially engage with the hard plastic.
Respironics	FitLife	S	Small (T101)	Excellent	The FitLife liner was developed specifically for the Respironics FitLife mask.
Respironics	FitLife	L/XL	Large/XL (T103)	Excellent	The FitLife liner was developed specifically for the Respironics FitLife mask.
Respironics	Full Life	S	Small (F101)	Poor	Quick Stick tabs do not reach hard plastic.
Respironics	Full Life	M	Medium (F102)	Poor	Quick Stick tabs do not reach hard plastic.
Respironics	Full Life	L	Large (F103)	Poor	Quick Stick tabs do not reach hard plastic.
Sleepnet	Ascend	S	Small (F101)	Excellent	Align the top edge of the liner with the top of the gel cushion.
Sleepnet	Ascend	M	Medium (F102)	Excellent	Align the top edge of the liner with the top of the gel cushion.
Sleepnet	Ascend	L	Large (F103)	Excellent	Align the top edge of the liner with the top of the gel cushion.
Sleepnet	Mojo	S	Small (F101)	Excellent	Align the top edge of the liner with the top of the gel cushion.
Sleepnet	Mojo	M	Medium (F102)	Excellent	Align the top edge of the liner with the top of the gel cushion.
Sleepnet	Mojo	L	Large (F103)	Excellent	Align the top edge of the liner with the top of the gel cushion.