

## The Best in Comfort and Technology



- Soft
- Light
- Quiet
- Easy to fit
- No skin irritation from plastic
- Comfortable seal on many face types
- Compatible with common CPAP/Bi-Level machines
- Durable
- Easy to clean
- Economical

**“Love the mask”**

*Raymond B.*

**“Excellent Mask,  
Very Comfortable”**

*Rich H.*

**“I fell in love with this  
one (SleepWeaver 3D)”**

*John B.*

**“Best I have used,  
works very well and  
comfortable.”**

*Mary Jane S.*



### **New wing-shaped stay:**

Routes the headgear under the ears, ensuring comfort and a great fit.

### **Permanent 360° Swivel:**

Provides reliable performance for the tubing in any position, including over the headboard.

### **Nothing Breathes Like Cloth:**

Breathable cloth design eliminates fogging common in hard plastic masks. Keeps skin dry and comfortable.

### **Lightweight:**

What could be lighter than a thin piece of cloth?

### **Accommodates Eye Glasses:**

For those who want to read or watch TV while in bed, SleepWeaver 3D allows users to wear eye glasses with the mask on.

### **The Joy of Soft Cloth:**

No one can deny preferring a soft cloth against their face as opposed to rubber, silicone and other gels. Cloth masks eliminate the need for a cushion or nasal pad.

### **A new and more flexible elastic interface material:**

This creates our best seal ever on many face types.

### **New rounded edge straps:**

Reduces strap marks.


# CPAPX



# SleepWeaver 3D ordering information

Product	Photos	HCPCS Codes	Part Number
3D Mask & Large Headgear		A7034-Mask A7035-Headgear	101491
3D Mask & Regular Headgear		A7034-Mask A7035-Headgear	101486
3D Mask Only		A7034-Mask	101496
3D Large Headgear Only		A7035-Headgear	101517
3D Regular Headgear Only		A7035-Headgear	101516
Feather Weight Tube Only		A7037-Tube	100477
Hose Clip		N/A	100775

## SleepWeaver 3D Sizing Guide



Watch SleepWeaver 3D Videos Online:  
[www.circadiancance.com/3dvideos](http://www.circadiancance.com/3dvideos)

101563 REV 1

Phone: 888-825-9640 or +1-724-858-2837  
[orders@circadiancance.com](mailto:orders@circadiancance.com) [www.circadiancance.com](http://www.circadiancance.com)



**SleepWeaver**  
 Soft Cloth PAP Masks