

## Deep restorative sleep for mouth breathers.

### The World's First Cloth Full Face PAP Mask



**CPAPXCHANGE**  
SLEEP WELL TONIGHT

#### **Super Lightweight:**

Users barely notice the weight of this mask.

#### **3 Sizes:**

Multiple size options allow users to find the perfect fit for their face; expands the number of users the mask will fit.

#### **Moisture Vapor Breathable:**

High-tech material allows the skin to breathe naturally, wicks away moisture and leaves no pressure marks common with rubber or silicone seals.

#### **Loose Fit:**

As with all SleepWeaver masks, the unique "loose fit" provides ultimate comfort and eliminates strap marks.

#### **Comfortable Seal for Mouth Breathers:**

The dual interface design individually seals around the nose and the mouth allowing only soft cloth to touch the skin.

#### **Open Face Headgear Design:**

Allows virtually unrestricted vision; great for people who like to fall asleep while reading or watching TV. Accommodates eyeglasses.

#### **90° Elbow Swivel Connector:**

Allows users to sleep in more natural positions; easy to remove and clean.

#### **Gentle Exhalation:**


Exhalation holes are in the front of the mask providing a quiet and gentle flow of air.

#### **Seals in a Snap:**












Two quick release snaps makes Anew one of the easiest masks to put on and take off.

#### **Easy to Clean:**

Quick & simple to clean the material that routinely comes in contact with the skin surface.



# SleepWeaver Anew Ordering Information

Product	Photos	Sizes	Part Numbers	Color Options
SleepWeaver Anew Mask & Headgear		Small Regular Large		
			Small 100955	Blue 
			Regular 100960	Blue 
			Large 100965	Blue 
SleepWeaver Anew Mask Only		Small Regular Large		
			Small 100970	Blue 
			Regular 100975	Blue 
			Large 100980	Blue 
SleepWeaver Anew Headgear Only		Small Regular		
			Small 101145	Black 
			Regular 101084	Black 

## SleepWeaver Anew Accessories

					
Tether Strap - 100638	Hose Clip - 100775	DVD 10 Pack - 101086	45 cm Tube - 100477	Washer - 100778	Connector - 100777
Tether Strap attaches to Tube & Headgear	Hose Clip prevents Tube movement	SleepWeaver Anew Instructional DVD video with printable Instructions for Use	SleepWeaver Feather Weight Tube	SleepWeaver Anew Washer Replacement Part	SleepWeaver Anew Threaded Connector Replacement Part



**EMERGO EUROPE**  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands  
+31 (0)70 345 8570

  
1060 Corporate Lane  
Export, PA 15632 USA  
USA: 888-825-9640  
+1-724-858-2837  
[www.Circadiance.com](http://www.Circadiance.com)



# SleepWeaver®

## The World's Only Line of Soft Cloth PAP Masks