

SleepWeaver® Anew™ Online Sizing Template

Instructions for Use:

1. To ensure the Sizing Guide is true to scale, verify the printer is set to print at 100%.
2. After printing, check the size of the printed version of the Sizing Guide by lining up a ruler with the ruler scale on the guide.
3. Carefully cut out the Sizing Guide on the dotted line.
4. Remember, a soft fabric mask is different than a hard plastic mask. Fabric can be pulled tighter and relaxes as it is used daily.
5. If you fall "in between" sizes, please order the smaller size of Anew.

