

# SperaChamber

**Anti-static valved inhaler chamber for metered dose inhaler**

SperaChamber Only

Item#  
**RES195**

SperaChamber with Dog Mask — Pediatric (3-6 Years)

Item#  
**RES195MS** Small  
**RES195MM** Medium



- Latex Free
- DEHP Free
- Universal adapter with spiral vents
- Mask made from soft and comfortable material which easily contours to the patient's face to minimize medication loss and reduce waste
- Mask's exhalation safety valve releases CO<sub>2</sub> when patient breathes out
- Chamber features one-way valve that opens to allow inhalation of aerosol and closes to prevent patient from exhaling into the chamber
- Drug delivery is further enhanced by the antistatic quality of the chamber, a feature that helps prevent medication from clinging to the chamber wall

## Specifications

**Chamber Volume** 160cm<sup>3</sup>

**Weight** 65 grams

**Mouthpiece Diameter** 22mm

**Dimensions** 53mm x 150mm

**Storage Temperature** 14°–113°F (-10°–45°C)

**Storage Relative Humidity** 30%–80%

**Single Patient Use** Yes



Adult size mask  
also available  
**RES195ML**

