



Meet the Transcend Micro™

Small but Mighty Powerhouse
is Packed with Features You'll Love



**Full-range,
full-featured
auto-CPAP**



**Smallest &
lightest CPAP on
the market**



**Whisper
quiet
operation**



**Compatible
with most
masks***



**Multiple
power options—
battery & solar**



**AirMist HME –
waterless
humidification**



**Built-in drying
cycle for moisture
removal**



**MySleepDash™
App therapy
tracking**



**AirRelief –
expiratory
pressure relief**



**GentleRise™
pressure
ramp feature**



**FAA compliant
for inflight use**



**2-year
warranty**



**Auto-adjusting
leak compensation**



**Bluetooth®
technology
enabled**

“The Transcend is very portable and convenient for travel. And just as comfortable as my regular-use machine. This machine packs easily in carry-on luggage.”

— Roger C.

* Compatible with any CPAP mask with a standard 22mm connection

Detailed Features and Benefits of the Transcend Micro

Portability: The small (3.6" x 3.6" x 2.4") size and light weight (0.48 lbs) make it easy to fit into a carry-on. It's ready to travel anywhere!

Sound Volume/Level: The whisper quiet design generates around 27 dBA, allowing for a more peaceful sleeping environment. Because we know how important sound is while using a CPAP device, we ship every device with a WhisperSoft™ muffler.

Comfort Settings: Featuring GentleRise™ technology designed to help you fall asleep more comfortably by delivering a lower pressure when you turn on the device, which gradually increases over time until the desired level is reached. Plus, it's adjustable, so you can fine-tune your ramp times or turn off the ramp feature to support your therapy needs. Additionally, the AirRelief feature is designed to help you breathe easier while you sleep by reducing pressure on exhalation.

Drying Mode: Removes excess moisture from your mask and tubing. Simply activate the 30-minute drying cycle upon waking to start a steady stream of air moving through the hose and mask.

Power Options: No matter where your adventures lead you, you'll be equipped with versatile power options.

- The small, robust Transcend P8 multi-night battery provides 16-20 hours of runtime in a compact, lightweight design.
- The durable, military-grade Transcend Portable Solar Battery Charger can recharge the P8 battery with help from the sun, even in low-light conditions.

Mask Compatibility: Have a favorite mask? No problem. The Transcend Micro can be used with any mask with a standard 22mm connection, eliminating the need for additional equipment.

Humidification: The optional Transcend AirMist™ HME (heat and moisture exchanger) provides pass over humidification without the need for a separate water chamber. It's all about comfort, without extra baggage.

Airline Travel: FAA-approved for in-flight use, allowing uninterrupted CPAP therapy during air travel.

App for Sleep Tracking: The MySleepDash™ app is like a CPAP's best friend. It provides information on therapy usage (usage hours, mask leak summary, events per hour monitoring and mask removal detection), comfort setting adjustments, 1-click troubleshooting and 1-click product feedback.

